


Implementing interventions and collecting data electronically is becoming the norm.



1. What capabilities do these technologies provide?
Which apply to your work?
2. Where are the trends moving?
3. How to design state-of-the-art/science interventions





High End
\$250 - \$3,500

Mid Range
\$61 - \$250

Low End
\$5 - \$60



Most MedSignals® features resulted from researcher needs



Programmable

- Online
- On the device
- **Via your laptop**
- Point and Click on the **Webpage.** The phone downloads new settings to the device from your desk.
- On the **device** using 3 buttons.
- From your own computer with a **phone simulator** (known as an I-Set), synchronize the device with your new settings

Main Menu | My Info | Bin Details | Device Settings | Pill-Taking History | Change Password | Sign Off

BIN PROGRAMMING

Patient Name: Kimberly Hahn
Patient Number: 7011001511
Patient Group No: 99999

Bin Number:

PILL DETAILS

1 Pill Name
Dosage (mg)
Number of Pills per Prescription
Prescription Fill Date

DEVICE DETAILS SETTINGS

2

Parameters	Current Setting (Last Upload: CST)	New Settings To Download
Effective Date	Enter New Settings Effective Date	<input type="text" value="5/16/2008"/>
Pill Cycle	Not Set	<input type="button" value="Bin Off"/>
Pill Type	Not Set	<input type="button" value="Half Pills"/>
Number Of Pills	Not Set	<input type="button" value="0.5"/>
Timing Format	Not Set	<input type="button" value="Regular Intervals"/>

PILL TIME SETTINGS

3

Current Dose Time(s)	New Dose Time(s)
Not Set	Not Applicable

WARNING SETTINGS

4

Current Warnings	New Warnings																																								
<table><thead><tr><th>Warning Description</th><th>Selected</th></tr></thead><tbody><tr><td>TAKE W/FOOD</td><td><input type="checkbox"/></td></tr><tr><td>TAKE W/MEALS</td><td><input type="checkbox"/></td></tr><tr><td>TAKE W/WATER</td><td><input type="checkbox"/></td></tr><tr><td>AVOID SUNLIGHT</td><td><input type="checkbox"/></td></tr><tr><td>EMPTY STOMACH</td><td><input type="checkbox"/></td></tr><tr><td>AVOID ALCOHOL</td><td><input type="checkbox"/></td></tr><tr><td>DROWSINESS</td><td><input type="checkbox"/></td></tr><tr><td>AVOID HI FAT</td><td><input type="checkbox"/></td></tr><tr><td>TAKE NO FOOD</td><td><input type="checkbox"/></td></tr></tbody></table>	Warning Description	Selected	TAKE W/FOOD	<input type="checkbox"/>	TAKE W/MEALS	<input type="checkbox"/>	TAKE W/WATER	<input type="checkbox"/>	AVOID SUNLIGHT	<input type="checkbox"/>	EMPTY STOMACH	<input type="checkbox"/>	AVOID ALCOHOL	<input type="checkbox"/>	DROWSINESS	<input type="checkbox"/>	AVOID HI FAT	<input type="checkbox"/>	TAKE NO FOOD	<input type="checkbox"/>	<table><thead><tr><th>Warning Description</th><th>Select</th></tr></thead><tbody><tr><td>TAKE W/FOOD</td><td><input type="checkbox"/></td></tr><tr><td>TAKE W/MEALS</td><td><input type="checkbox"/></td></tr><tr><td>TAKE W/WATER</td><td><input type="checkbox"/></td></tr><tr><td>AVOID SUNLIGHT</td><td><input type="checkbox"/></td></tr><tr><td>EMPTY STOMACH</td><td><input type="checkbox"/></td></tr><tr><td>AVOID ALCOHOL</td><td><input type="checkbox"/></td></tr><tr><td>DROWSINESS</td><td><input type="checkbox"/></td></tr><tr><td>AVOID HI FAT</td><td><input type="checkbox"/></td></tr><tr><td>TAKE NO FOOD</td><td><input type="checkbox"/></td></tr></tbody></table>	Warning Description	Select	TAKE W/FOOD	<input type="checkbox"/>	TAKE W/MEALS	<input type="checkbox"/>	TAKE W/WATER	<input type="checkbox"/>	AVOID SUNLIGHT	<input type="checkbox"/>	EMPTY STOMACH	<input type="checkbox"/>	AVOID ALCOHOL	<input type="checkbox"/>	DROWSINESS	<input type="checkbox"/>	AVOID HI FAT	<input type="checkbox"/>	TAKE NO FOOD	<input type="checkbox"/>
Warning Description	Selected																																								
TAKE W/FOOD	<input type="checkbox"/>																																								
TAKE W/MEALS	<input type="checkbox"/>																																								
TAKE W/WATER	<input type="checkbox"/>																																								
AVOID SUNLIGHT	<input type="checkbox"/>																																								
EMPTY STOMACH	<input type="checkbox"/>																																								
AVOID ALCOHOL	<input type="checkbox"/>																																								
DROWSINESS	<input type="checkbox"/>																																								
AVOID HI FAT	<input type="checkbox"/>																																								
TAKE NO FOOD	<input type="checkbox"/>																																								
Warning Description	Select																																								
TAKE W/FOOD	<input type="checkbox"/>																																								
TAKE W/MEALS	<input type="checkbox"/>																																								
TAKE W/WATER	<input type="checkbox"/>																																								
AVOID SUNLIGHT	<input type="checkbox"/>																																								
EMPTY STOMACH	<input type="checkbox"/>																																								
AVOID ALCOHOL	<input type="checkbox"/>																																								
DROWSINESS	<input type="checkbox"/>																																								
AVOID HI FAT	<input type="checkbox"/>																																								
TAKE NO FOOD	<input type="checkbox"/>																																								

Notify, if Overdose/Underdose for consecutive day(s) ☒




When an analog phone line isn't available in your office, the **MedSignals® I-Set™** provides quick office set-ups.



- Bypass phone, uploads straight to server when computer is online.
- Program many devices for a protocol in minutes.
- Download devices in office visits, without waiting for phone dial-ups.
- Access website globally without int'l phone costs.



When patients need monitoring on more than 4 drugs, multi-packs can be integrated into common data reports.




Trends in Home-Based Monitoring:

- Instant alerts, notifications to care team
- Wireless (embedded cell phone chips) in devices – not Bluetooth-ing to individuals' cell phones.
- Co-morbidity drugs with reports to multiple


nubs (MedSignals, AI & I Actuarious, etc.)

- Nurse-led Call Centers for exception handling, barriers management



Opportunities for Researchers:

- **Co-morbidities: comprehensive disease management**
- **Non-standard use of bins: methadone weaning**
- **Event reminders beyond meds: “virtual” bins**
- **Non-pill medication forms (blister packs, injections, refrigerated drugs)**
- **Solar-powered devices (homeless, African HIV)**
- **Cost-effectiveness and comparative effectiveness measures**
- **Non-English languages**
- **Call center and phone IVR feedback loop**
- **Apps for smart phones**
- **RFID-tagged drugs and other diversion technologies**





Contact Information:

Vesta Brue

Founder/CEO

vesta@medsignals.com

Phone: 210-222-2067

MedSignals LLC

217 Alamo Plaza

San Antonio, TX 78205




MedSignals® Functions in Intervention Mode

- Records each pill-taking event in memory
- Signals at pill times in 4 multi-sensory ways (**beeps, flashing, text, VOICE**).
- Voice announces how many pills to take and how to take them (9 optional warnings) in **English, Spanish** (others coming.)
- Provides feedback at the press of a button on today's usage (last, next, and how many pills today) in the Intervention mode.
- Uploads lid opening data to server database and imports to your database **or the patient's Electronic Medical Record**
- Displays Smart Charts on patient's secured, customized web page
- Dashboards and aggregate charts help study managers **quickly respond to non-adherents** (using your own criteria) for exception handling.



MedSignals® Functions in Control Condition Mode


- Turn off all signals and patient feedback at the push of a button.
Tracks natural pill-taking patterns thereafter.
- Records each lid-opening event in memory
- Automatically uploads lid opening data to server database, imports to your database or the patient's electronic record
- Charts activity on participant's secured web page to quickly spot malfunctions or participant attrition. Dashboards and aggregate charts help study managers quickly respond to non-adherents (using your own criteria).



Unique Benefits to Researchers:

- Improve data collection reliability and ease. Reduce attrition. Identify and remedy non-performing participants early. Export data anytime, any way. .
- Conduct polydrug studies efficiently. Track multiple drugs simultaneously.
- Set devices remotely. Two-way communications lets you change settings when prescriptions change, without an office visit.
- Reuse devices with new participants, new studies. Rechargeable battery should last years.

>>> more >>>>



- Convert quickly from Intervention to control condition data collection
- Bluetooth connectivity from other home-based bio-monitors, capturing real-time metrics.
- Cut costs, labor, and participant burden.
- Continuous data collection enables adaptive trials.
- Execute customized study interventions with reminder signals.
- Match device models and functions to your study protocols.
- Choose frequency of uploads to match your needs and budget.

